



Work Instruction for: CGCR Cleaning Protocols

Note: Sanitizer used for this protocol is a Alkyl dimethyl benzyl ammonium chloride cleaner. Brand names include Barbicide and Lobacten. Dilute following the directions on the label.

A. Cages /carriers

1. Remove all bedding, dishes, and other items. Sweep out all debris, clean out any feces or droppings out as well. All organic matter (poo, pee, food) be bleach wiped /scrubbed clean and removed before the sterilization process.
2. Once the cage is cleared out, spray with Sanitizer in pump sprayer all surfaces the inside of cages or kennel and the doors both inside and out. Allow to dry completely for 1/2 hour.
3. All surfaces must be bleach wiped out with warm bleach water – use 3 cups bleach per gallon of water.

B. Floors, Walls and Other Surfaces

1. Sweep up all debris, clean out any feces or droppings out as well. Use bleach water to remove all organic matter.
2. Spray with Sanitizer. Allow to dry completely for 1/2 hour.
3. All surfaces must be bleach wiped out with warm bleach water (dilution of 2 cups bleach per gallon of water).

C. Litter Scoops, Scoop Holders, Dishes, Washable Toys

1. Remove as much organic matter by rinsing and wiping.
2. Run through dishwasher using sterilizing mode. Small litter boxes are best done in dishwasher as well.
3. Large litter boxes can be cleaned by hand using the same method as carriers. Very important to include all sides, interior and exterior and under lip. Boil water in a kettle to help loosen and dry-stuck feces.

Note: non-washable/sterilizable toys and cardboard scratchers must either go with the cats or be thrown away.